



## NF-1 Polyclonal Antibody

Cat No: HR1AP7478

For research use only

### Overview

Product Name	NF-1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GenID?Human?	NFIA
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NFIA; KIAA1439; Nuclear factor 1 A-type; NF1-A; Nuclear factor 1A; CCAAT-box-binding transcription factor; CTF; Nuclear factor I/A; NF-I/A; NFI-A; TGGCA-binding protein; NFIB; Nuclear factor 1 B-type
Background	nuclear factor I A(NFIA) Homo sapiens This gene encodes a member of the NF1 (nuclear factor 1) family of transcription factors. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2011],