



# OSR1 (phospho Thr185) Polyclonal Antibody

Cat No: HR1AP9251

For research use only

## Overview

Product Name	OSR1 (phospho Thr185) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	65kDa
GenID?Human?	OXSR1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	OXSR1; KIAA1101; OSR1; Serine/threonine-protein kinase OSR1; Oxidative stress-responsive 1 protein
Background	oxidative stress responsive 1(OXSR1) Homo sapiens The product of this gene belongs to the Ser/Thr protein kinase family of proteins. It regulates downstream kinases in response to environmental stress, and may play a role in regulating the actin cytoskeleton. [provided by RefSeq, Jul 2008],