



FoxL2 (phospho Ser263) Polyclonal Antibody

Cat No: HR1AP8372

For research use only

Overview

Product Name	FoxL2 (phospho Ser263) Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40kDa
GeneID?Human?	FOXL2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FOXL2; Forkhead box protein L2
Background	forkhead box L2(FOXL2) Homo sapiens This gene encodes a forkhead transcription factor. The protein contains a fork-head DNA-binding domain and may play a role in ovarian development and function. Expansion of a polyalanine repeat region and other mutations in this gene are a cause of blepharophimosis syndrome and premature ovarian failure 3. [provided by RefSeq, Jul 2016],