

A20 Polyclonal Antibody

Cat No: HR1AP9474

For research use only

Overview

Product Name	A20 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	95kDa
GeneID?Human?	ANKRD20A1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ANKRD20A1; ANKRD20A; Ankyrin repeat domain-containing protein 20A1
Background	similarity:Contains 5 ANK repeats.,tissue specificity:Widely expressed; not found in kidney, nor in liver.,