



## TEF-4 Polyclonal Antibody

Cat No: HR1AP8878

For research use only

### Overview

Product Name	TEF-4 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GenID?Human?	TEAD2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TEAD2; TEF4; Transcriptional enhancer factor TEF-4; TEA domain family member 2; TEAD-2
Background	function:Putative transcription factor that binds to the SPH and GT-IIC "enhancers" (5'-GTGGAATGT-3'). May be involved in the gene regulation of neural development. Binds to the M-CAT motif.,similarity:Contains 1 TEA DNA-binding domain.,