



## PHF3 Polyclonal Antibody

Cat No: HR1AP4303

For research use only

### Overview

Product Name	PHF3 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	230kDa
GenID?Human?	PHF3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PHF3; KIAA0244; PHD finger protein 3
Background	PHD finger protein 3(PHF3) Homo sapiens This gene encodes a member of a PHD finger-containing gene family. This gene may function as a transcription factor and may be involved in glioblastomas development. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2014],