



# PHF3 Polyclonal Antibody

Cat No: HR1AP4303

For research use only

## Overview

|                       |                                                                                                                                                                                                                                                                                                             |
|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | PHF3 Polyclonal Antibody                                                                                                                                                                                                                                                                                    |
| Source                | Rabbit                                                                                                                                                                                                                                                                                                      |
| Applications          | WB,IHC-p,ELISA                                                                                                                                                                                                                                                                                              |
| Species Reactivity    | Human,Mouse                                                                                                                                                                                                                                                                                                 |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                                             |
| Immunogen             |                                                                                                                                                                                                                                                                                                             |
| Species               | Rabbit                                                                                                                                                                                                                                                                                                      |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                                                |
| Isotype               |                                                                                                                                                                                                                                                                                                             |
| Clonality             |                                                                                                                                                                                                                                                                                                             |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                                                     |
| Observed band         | 230kDa                                                                                                                                                                                                                                                                                                      |
| GeneID?Human?         | PHF3                                                                                                                                                                                                                                                                                                        |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                                             |
| Cellular localization |                                                                                                                                                                                                                                                                                                             |
| Alternative Names     | PHF3; KIAA0244; PHD finger protein 3                                                                                                                                                                                                                                                                        |
| Background            | PHD finger protein 3(PHF3) Homo sapiens This gene encodes a member of a PHD finger-containing gene family. This gene may function as a transcription factor and may be involved in glioblastomas development. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2014], |