

Fos B Polyclonal Antibody

Cat No: HR1AP3470

For research use only

Overview

Product Name	Fos B Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	35kDa
GeneID?Human?	FOSB
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FOSB; G0S3; Protein fosB; G0/G1 switch regulatory protein 3
Background	FosB proto-oncogene, AP-1 transcription factor subunit(FOSB) Homo sapiens The Fos gene family consists of 4 members: FOS, FOSB, FOSL1, and FOSL2. These genes encode leucine zipper proteins that can dimerize with proteins of the JUN family, thereby forming the transcription factor complex AP-1. As such, the FOS proteins have been implicated as regulators of cell proliferation, differentiation, and transformation. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008],