

Aquaporin 4 Mouse mAb

Cat No: HR1AM2173

For research use only

Overview

Product Name	Aquaporin 4 Mouse mAb
Source	Mouse
Applications	WB, IHC, IF
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000 IF 1:100-200 IHC 1:200-500
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	48kDa
GeneID?Human?	361
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	AQP4, aquaporin 4, NMO-IgG
Background	Aquaporin 4 belongs to the Aquaporin family of integral membrane proteins that conduct water through the cell membrane.