



## CBX6 Polyclonal Antibody

Cat No: HR1AP2995

For research use only

### Overview

Product Name	CBX6 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	44kDa
GenID?Human?	CBX6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CBX6; Chromobox protein homolog 6
Background	function:Involved in maintaining the transcriptionally repressive state of genes. Modifies chromatin, rendering it heritably changed in its expressibility.,similarity:Contains 1 chromo domain.,subunit:Component of the chromatin-associated polycomb complex (PcG),,