



AF-10 Polyclonal Antibody

Cat No: HR1AP2720

For research use only

Overview

Product Name	AF-10 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	109kDa
GeneID?Human?	MLLT10
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MLLT10; AF10; Protein AF-10; ALL1-fused gene from chromosome 10 protein
Background	myeloid/lymphoid or mixed-lineage leukemia; translocated to, 10(MLLT10) Homo sapiens This gene encodes a transcription factor and has been identified as a partner gene involved in several chromosomal rearrangements resulting in various leukemias. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2010],