

Sp1 (phospho Thr739) Polyclonal Antibody

Cat No: HR1AP2514

For research use only

Overview

Product Name	Sp1 (phospho Thr739) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	81kDa
GeneID?Human?	SP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SP1; TSFP1; Transcription factor Sp1
Background	Sp1 transcription factor(SP1) Homo sapiens The protein encoded by this gene is a zinc finger transcription factor that binds to GC-rich motifs of many promoters. The encoded protein is involved in many cellular processes, including cell differentiation, cell growth, apoptosis, immune responses, response to DNA damage, and chromatin remodeling. Post-translational modifications such as phosphorylation, acetylation, glycosylation, and proteolytic processing significantly affect the activity of this protein, which can be an activator or a repressor. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2014],