



GAL4 Activation Domain Polyclonal Antibody

Cat No: HR1AP9341

For research use only

Overview

Product Name	GAL4 Activation Domain Polyclonal Antibody
Source	Rabbit
Applications	WB
Species Reactivity	
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	
Observed band	kDa
GenID?Human?	
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>The GAL4 protein of <i>Saccharomyces cerevisiae</i> is one of the most thoroughly characterized transcriptional activators. Since the N-terminal 147 amino acid residues of GAL4 are sufficient to mediate specific and strong binding to DNA, but are incapable of efficient transcriptional activation, this protein fragment has frequently been used to confer specific DNA binding in experiments examining transcriptional activation functions of heterologous proteins. This approach is facilitated by the finding that higher eukaryotes lack endogenous proteins that enhance transcription from the consensus GAL4-binding site. Fusions between GAL4 (aa 1-147) and activating domains from a variety of transcriptional regulatory proteins can activate transcription in yeast, plant, insects and mammalian cells. A unique "two-hybrid" system has been developed using GAL4 fusions in yeast to identify specific protein-protein interactions.</p>