

CDYL2 Polyclonal Antibody

Cat No: HR1AP3058

For research use only

Overview

Product Name	CDYL2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ICC,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	56kDa
GeneID?Human?	CDYL2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CDYL2; Chromodomain Y-like protein 2; CDY-like 2
Background	similarity:Contains 1 chromo domain.,tissue specificity:Ubiquitously expressed.,