

HMG-1 (Acetyl Lys12) Polyclonal Antibody

Cat No: HR1AP2246

For research use only

Overview

Product Name	HMG-1 (Acetyl Lys12) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	25kDa
GeneID?Human?	HMGB1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	HMGB1; HMG1; High mobility group protein B1; High mobility group protein 1; HMG-1
Background	high mobility group box 1(HMGB1) Homo sapiens This gene encodes a protein that belongs to the High Mobility Group-box superfamily. The encoded non-histone, nuclear DNA-binding protein regulates transcription, and is involved in organization of DNA. This protein plays a role in several cellular processes, including inflammation, cell differentiation and tumor cell migration. Multiple pseudogenes of this gene have been identified. Alternative splicing results in multiple transcript variants that encode the same protein. [provided by RefSeq, Sep 2015],