



TCF-3 Polyclonal Antibody

Cat No: HR1AP8793

For research use only

Overview

Product Name	TCF-3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	63kDa
GenelD?Human?	TCF7L1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TCF7L1; TCF3; Transcription factor 7-like 1; HMG box transcription factor 3; TCF-3
Background	transcription factor 7 like 1(TCF7L1) <i>Homo sapiens</i> This gene encodes a member of the T cell factor/lymphoid enhancer factor family of transcription factors. These transcription factors are activated by beta catenin, mediate the Wnt signaling pathway and are antagonized by the transforming growth factor beta signaling pathway. The encoded protein contains a high mobility group-box DNA binding domain and participates in the regulation of cell cycle genes and cellular senescence. [provided by RefSeq, Nov 2010],