



## ERF (phospho Thr526) Polyclonal Antibody

Cat No: HR1AP6279

For research use only

### Overview

Product Name	ERF (phospho Thr526) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	58kDa
GenelD?Human?	ERF
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ERF; ETS domain-containing transcription factor ERF; Ets2 repressor factor; PE-2
Background	ETS2 repressor factor(ERF) Homo sapiens ETS2 is a transcription factor and protooncogene involved in development, apoptosis, and the regulation of telomerase. The protein encoded by this gene binds to the ETS2 promoter and is a strong repressor of ETS2 transcription. Several transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Aug 2015],