



SCML2 Polyclonal Antibody

Cat No: HR1AP4942

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | SCML2 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 77kDa |
| GenID?Human? | SCML2 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | SCML2; Sex comb on midleg-like protein 2 |
| Background | sex comb on midleg-like 2 (Drosophila)(SCML2) Homo sapiens This gene encodes a member of the Polycomb group proteins. These proteins form the Polycomb repressive complexes which are involved in transcriptional repression. The encoded protein binds histone peptides that are monomethylated at lysine residues and may be involved in regulating homeotic gene expression during development. [provided by RefSeq, Jun 2010], |