

SCML2 Polyclonal Antibody

Cat No: HR1AP4942

For research use only

Overview

Product Name	SCML2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	77kDa
GeneID?Human?	SCML2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SCML2; Sex comb on midleg-like protein 2
Background	sex comb on midleg-like 2 (Drosophila)(SCML2) Homo sapiens This gene encodes a member of the Polycomb group proteins. These proteins form the Polycomb repressive complexes which are involved in transcriptional repression. The encoded protein binds histone peptides that are monomethylated at lysine residues and may be involved in regulating homeotic gene expression during development. [provided by RefSeq, Jun 2010],