

EF-1 ? Polyclonal Antibody

Cat No: HR1AP3340

For research use only

Overview

Product Name	EF-1 ? Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GeneID?Human?	EEF1G
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	EEF1G; EF1G; Elongation factor 1-gamma; EF-1-gamma; eEF-1B gamma
Background	eukaryotic translation elongation factor 1 gamma(EEF1G) Homo sapiens This gene encodes a subunit of the elongation factor-1 complex, which is responsible for the enzymatic delivery of aminoacyl tRNAs to the ribosome. This subunit contains an N-terminal glutathione transferase domain, which may be involved in regulating the assembly of multisubunit complexes containing this elongation factor and aminoacyl-tRNA synthetases. [provided by RefSeq, Jul 2008],