



## RRS1 Polyclonal Antibody

Cat No: HR1AP4505

For research use only

### Overview

|                       |   |
|-----------------------|---|
| Product Name          | RRS1 Polyclonal Antibody  |
| Source                | Rabbit  |
| Applications          | WB,IHC-p,ELISA  |
| Species Reactivity    | Human,Mouse,Rat   |
| Recommended Dilutions |   |
| Immunogen             |   |
| Species               | Rabbit  |
| Storage               | -20°C/1 year  |
| Isotype               |   |
| Clonality             |   |
| Concentration         | 1 mg/ml   |
| Observed band         | 41kDa   |
| GenID?Human?          | RRS1  |
| Human Swiss-Prot No.  |   |
| Cellular localization |   |
| Alternative Names     | RRS1; KIAA0112; RRR; Ribosome biogenesis regulatory protein homolog               |
| Background            | function:Involved in ribosome biogenesis.,similarity:Belongs to the RRS1 family., |