

RRS1 Polyclonal Antibody

Cat No: HR1AP4505

For research use only

Overview

Product Name	RRS1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	41kDa
GeneID?Human?	RRS1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RRS1; KIAA0112; RRR; Ribosome biogenesis regulatory protein homolog
Background	function:Involved in ribosome biogenesis.,similarity:Belongs to the RRS1 family.,