



CCDC109A Polyclonal Antibody

Cat No: HR1AP9055

For research use only

Overview

Product Name	CCDC109A Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	MCU
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MCU; C10orf42; CCDC109A; Calcium uniporter protein; mitochondrial; Coiled-coil domain-containing protein 109A
Background	mitochondrial calcium uniporter(MCU) Homo sapiens This gene encodes a calcium transporter that localizes to the mitochondrial inner membrane. The encoded protein interacts with mitochondrial calcium uptake 1. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2012].