

IRF-3 (phospho Ser385) Polyclonal Antibody

Cat No: HR1AP7034

For research use only

Overview

Product Name	IRF-3 (phospho Ser385) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	57kDa
GeneID?Human?	IRF3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	IRF3; Interferon regulatory factor 3; IRF-3
Background	interferon regulatory factor 3(IRF3) Homo sapiens This gene encodes a member of the interferon regulatory transcription factor (IRF) family. The encoded protein is found in an inactive cytoplasmic form that upon serine/threonine phosphorylation forms a complex with CREBBP. This complex translocates to the nucleus and activates the transcription of interferons alpha and beta, as well as other interferon-induced genes. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Nov 2011],