



## FUNDC1 Polyclonal Antibody

Cat No: HR1AP5494

For research use only

### Overview

Product Name	FUNDC1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	17kDa
GenID?Human?	FUNDC1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FUNDC1; FUN14 domain-containing protein 1
Background	FUN14 domain containing 1(FUNDC1) Homo sapiens This gene encodes a protein with a FUN14 superfamily domain. The function of the encoded protein is not known. [provided by RefSeq, Sep 2011],

