



NF-?B p65 Polyclonal Antibody

Cat No: HR1AP9599

For research use only

Overview

Product Name	NF-?B p65 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	60kDa
GenID?Human?	RELA NFKB3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	v-rel reticuloendotheliosis viral oncogene homolog A (avian)
Background	RELA proto-oncogene, NF-?B subunit(RELA) Homo sapiens NF-kappa-B is a ubiquitous transcription factor involved in several biological processes. It is held in the cytoplasm in an inactive state by specific inhibitors. Upon degradation of the inhibitor, NF-kappa-B moves to the nucleus and activates transcription of specific genes. NF-kappa-B is composed of NFKB1 or NFKB2 bound to either REL, RELA, or RELB. The most abundant form of NF-kappa-B is NFKB1 complexed with the product of this gene, RELA. Four transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2011].