

?-actin(HRP Conjugated) Mouse mAb Cat No: HR1AM2149

For research use only

Overview

Product Name	?-actin(HRP Conjugated) Mouse mAb
Source	Mouse
Applications	WB
Species Reactivity	Human,Rat,Mouse,Monkey,Dog,Chicken,Hamster,Rabbit,Insect
Recommended Dilutions	WB 1:2,000-5,000 IHC 1:100
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	43kDakDa
GeneID?Human?	60
Human Swiss- Prot No.	
Cellular localization	
Alternative Names	actin antibody, b-actin antibody, beta-actin antibody
Background	?-actin is one of six different actin isoforms that have been identified. The actin molecules found in cells of various species and tissues tend to be very similar in their immunological and physical properties. Therefore, antibodies against ?-actin are useful as loading controls for Western Blotting. However it should be noted that levels of ?-actin may not be stable in certain cells. For example, expression of ?-actin in adipose tissue is very low and therefore ?-actin should not be used as loading control for these tissues.