



## BCoR Polyclonal Antibody

Cat No: HR1AP2874

For research use only

### Overview

Product Name	BCoR Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	192kDa
GenID?Human?	BCOR
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	BCOR; KIAA1575; BCL-6 corepressor; BCoR
Background	BCL6 corepressor(BCOR) Homo sapiens The protein encoded by this gene was identified as an interacting corepressor of BCL6, a POZ/zinc finger transcription repressor that is required for germinal center formation and may influence apoptosis. This protein selectively interacts with the POZ domain of BCL6, but not with eight other POZ proteins. Specific class I and II histone deacetylases (HDACs) have been shown to interact with this protein, which suggests a possible link between the two classes of HDACs. Several transcript variants encoding different isoforms have been found for this gene. A pseudogene of this gene is found on chromosome Y.[provided by RefSeq, Jun 2010],