



## YEATS2 Polyclonal Antibody

Cat No: HR1AP4831

For research use only

### Overview

Product Name	YEATS2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	150kDa
GenelD?Human?	YEATS2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	YEATS2; KIAA1197; YEATS domain-containing protein 2
Background	YEATS domain containing 2(YEATS2) Homo sapiens YEATS2 is a scaffolding subunit of the ADA2A (TADA2A; MIM 602276)-containing (ATAC) histone acetyltransferase complex (Wang et al., 2008 [PubMed 18838386]).[supplied by OMIM, Apr 2010],

