



## FOXP3 Polyclonal Antibody

Cat No: HR1AP5287

For research use only

### Overview

Product Name	FOXP3 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	47kDa
GenelD?Human?	FOXP3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FOXP3; IPEX; JM2; Forkhead box protein P3; Scurfin
Background	forkhead box P3(FOXP3) Homo sapiens The protein encoded by this gene is a member of the forkhead/winged-helix family of transcriptional regulators. Defects in this gene are the cause of immunodeficiency polyendocrinopathy, enteropathy, X-linked syndrome (IPEX), also known as X-linked autoimmunity-immunodeficiency syndrome. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008],