

Histone H1 (Di Methyl Lys25) Polyclonal Antibody

Cat No: HR1AP2327

For research use only

Overview

Product Name	Histone H1 (Di Methyl Lys25) Polyclonal Antibody
Source	Rabbit
Applications	WB
Species Reactivity	Human, Mouse, Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	
Observed band	17-25kDa
GeneID? Human?	HIST1H1A/HIST1H1C/HIST1H1D/HIST1H1E
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	H1K25ME2; HIST1H1A; H1F1; Histone H1.1; Histone H1a; HIST1H1C; H1F2; Histone H1.2; Histone H1c; Histone H1d; Histone H1s-1; HIST1H1D; H1F3; Histone H1.3; Histone H1c; Histone H1s-2; HIST1H1E; H1F4; Hi
Background	histone cluster 1 H1 family member a(HIST1H1A) Homo sapiens Histones are basic nuclear proteins responsible for nucleosome structure of the chromosomal fiber in eukaryotes. Two molecules of each of the four core histones (H2A, H2B, H3, and H4) form an octamer, around which approximately 146 bp of DNA is wrapped in repeating units, called nucleosomes. The linker histone, H1, interacts with linker DNA between nucleosomes and functions in the compaction of chromatin into higher order structures. This gene is intronless and encodes a replication-dependent histone that is a member of the histone H1 family. Transcripts from this gene lack polyA tails but instead contain a palindromic termination element. This gene is found in the large histone gene cluster on chromosome 6. [provided by RefSeq, Aug 2015],