



Strep-Tag Mouse mAb

Cat No: HR1AM2126

For research use only

Overview

Product Name	Strep-Tag Mouse mAb
Source	Mouse
Applications	WB
Species Reactivity	N/A
Recommended Dilutions	WB 1:5,000
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	N/AkDa
GeneID?Human?	N/A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	N/A
Background	The Strep-tag is a synthetic peptide consisting of eight amino acids: NWSHPQFEK (Strep tag) with high specificity and affinity towards streptavidin. The Strep tag antibody is a useful tool in analysis and affinity purification of Strep tag fusion proteins.