



USP32 Polyclonal Antibody

Cat No: HR1AP8884

For research use only

Overview

Product Name	USP32 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p, WB, ELISA
Species Reactivity	Human, Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	180kDa
GeneID?Human?	USP32
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	USP32; USP10; Ubiquitin carboxyl-terminal hydrolase 32; Deubiquitinating enzyme 32; Renal carcinoma antigen NY-REN-60; Ubiquitin thioesterase 32; Ubiquitin-specific-processing protease 32
Background	catalytic activity: Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol., similarity: Belongs to the peptidase C19 family., similarity: Contains 1 DUSP domain., similarity: Contains 3 EF-hand domains.,