




UBTD1 Polyclonal Antibody

Cat No: HR1AP8718

For research use only

Overview

Product Name	UBTD1 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	26kDa
GenID?Human?	UBTD1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	UBTD1; Ubiquitin domain-containing protein 1
Background	ubiquitin domain containing 1(UBTD1) Homo sapiens The degradation of many proteins is carried out by the ubiquitin pathway in which proteins are targeted for degradation by covalent conjugation of the polypeptide ubiquitin. This gene encodes a protein that belongs to the ubiquitin family of proteins. The encoded protein is thought to regulate E2 ubiquitin conjugating enzymes belonging to the UBE2D family. [provided by RefSeq, Mar 2014],



bioelsa