



PLC-XD1 Polyclonal Antibody

Cat No: HR1AP7834

For research use only

Overview

Product Name	PLC-XD1 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	PLCXD1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	PLCXD1; PI-PLC X domain-containing protein 1
Background	phosphatidylinositol specific phospholipase C X domain containing 1(PLCXD1) Homo sapiens This gene is the most terminal protein-coding gene in the pseudoautosomal (PAR) region on chromosomes X and Y. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Mar 2010],