

## T7-Tag Mouse mAb Cat No: HR1AM2122

For research use only

## **Overview**

Product Name	T7-Tag Mouse mAb
Source	Mouse
Applications	WB
Species Reactivity	N/A
Recommended Dilutions	WB 1:5,000
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	N/AkDa
GeneID?Human?	N/A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	N/A
Background	The T7 tag serves as a tag in many expression vectors including the pET system that is based on the very efficient T7 RNA polymerase expression system. The T7 tag is an 11 amino acid peptide encoded in the leader sequence of T7 bacteriophage gene10. This gene encodes a T7 major capsid protein.