



## Oxr1 Polyclonal Antibody

Cat No: HR1AP7815

For research use only

### Overview

|                       |   |
|-----------------------|---|
| Product Name          | Oxr1 Polyclonal Antibody  |
| Source                | Rabbit  |
| Applications          | WB,IHC-p,ELISA  |
| Species Reactivity    | Human,Mouse,Rat   |
| Recommended Dilutions |   |
| Immunogen             |   |
| Species               | Rabbit  |
| Storage               | -20°C/1 year  |
| Isotype               |   |
| Clonality             |   |
| Concentration         | 1 mg/ml   |
| Observed band         | 85kDa   |
| GenID?Human?          | OXR1  |
| Human Swiss-Prot No.  |   |
| Cellular localization |   |
| Alternative Names     | OXR1; Nbla00307; Oxidation resistance protein 1   |
| Background            | function:May be involved in protection from oxidative damage.,induction:By heat and oxidative stress.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the OXR1 family.,similarity:Contains 1 GRAM domain.,similarity:Contains 1 LysM repeat.,similarity:Contains 1 TLD domain., |