



NDUFC1 Polyclonal Antibody

Cat No: HR1AP7456

For research use only

Overview

Product Name	NDUFC1 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GenelID?Human?	NDUFC1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NDUFC1; NADH dehydrogenase [ubiquinone] 1 subunit C1; mitochondrial; Complex I-KFY1; CI-KFY1; NADH-ubiquinone oxidoreductase KFY1 subunit
Background	NADH:ubiquinone oxidoreductase subunit C1(NDUFC1) Homo sapiens The encoded protein is a subunit of the NADH:ubiquinone oxidoreductase (complex I), the first enzyme complex in the electron transport chain located in the inner mitochondrial membrane. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2010],