

# MARCH4 Polyclonal Antibody

**Cat No: HR1AP8032**

For research use only

## Overview

|                       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | MARCH4 Polyclonal Antibody                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Source                | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Applications          | IHC-p,ELISA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Species Reactivity    | Human,Mouse                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Immunogen             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Species               | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Isotype               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Clonality             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Observed band         | kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| GeneID?Human?         | Mar-04                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Cellular localization |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Alternative Names     | MARCH4; KIAA1399; RNF174; E3 ubiquitin-protein ligase MARCH4; Membrane-associated RING finger protein 4; Membrane-associated RING-CH protein IV; MARCH-IV; RING finger protein 174                                                                                                                                                                                                                                                                                                                                      |
| Background            | membrane associated ring-CH-type finger 4(MARCH4) Homo sapiens MARCH4 is a member of the MARCH family of membrane-bound E3 ubiquitin ligases (EC 6.3.2.19). MARCH enzymes add ubiquitin (see MIM 191339) to target lysines in substrate proteins, thereby signaling their vesicular transport between membrane compartments. MARCH4 reduces surface accumulation of several membrane glycoproteins by directing them to the endosomal compartment (Bartee et al., 2004 [PubMed 14722266]).[supplied by OMIM, Apr 2010], |