

HIBADH Polyclonal Antibody

Cat No: HR1AP5711

For research use only

Overview

Product Name	HIBADH Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	HIBADH
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	HIBADH; 3-hydroxyisobutyrate dehydrogenase; mitochondrial; HIBADH
Background	3-hydroxyisobutyrate dehydrogenase(HIBADH) Homo sapiens This gene encodes a mitochondrial 3-hydroxyisobutyrate dehydrogenase enzyme. The encoded protein plays a critical role in the catabolism of L-valine by catalyzing the oxidation of 3-hydroxyisobutyrate to methylmalonate semialdehyde. [provided by RefSeq, Nov 2011],