

GNPAT Polyclonal Antibody

Cat No: HR1AP8869

For research use only

Overview

Product Name	GNPAT Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	77kDa
GeneID?Human?	GNPAT
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GNPAT; DAPAT; DHAPAT; Dihydroxyacetone phosphate acyltransferase; DAP-AT; DHAP-AT; Acyl-CoA:dihydroxyacetonephosphateacyltransferase; Glycerone-phosphate O-acyltransferase
Background	glyceronephosphate O-acyltransferase(GNPAT) Homo sapiens This gene encodes an enzyme located in the peroxisomal membrane which is essential to the synthesis of ether phospholipids. Mutations in this gene are associated with rhizomelic chondrodysplasia punctata. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2015],