

# DHRS1 Polyclonal Antibody

**Cat No:** HR1AP5769

For research use only

## Overview

|                       |                                                                                                                                                                                                                                                                                                                                                                   |
|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | DHRS1 Polyclonal Antibody                                                                                                                                                                                                                                                                                                                                         |
| Source                | Rabbit                                                                                                                                                                                                                                                                                                                                                            |
| Applications          | IHC-p,ELISA                                                                                                                                                                                                                                                                                                                                                       |
| Species Reactivity    | Human,Mouse,Rat                                                                                                                                                                                                                                                                                                                                                   |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                                                                                                   |
| Immunogen             |                                                                                                                                                                                                                                                                                                                                                                   |
| Species               | Rabbit                                                                                                                                                                                                                                                                                                                                                            |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                      |
| Isotype               |                                                                                                                                                                                                                                                                                                                                                                   |
| Clonality             |                                                                                                                                                                                                                                                                                                                                                                   |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                           |
| Observed band         | kDa                                                                                                                                                                                                                                                                                                                                                               |
| GeneID?Human?         | DHRS1                                                                                                                                                                                                                                                                                                                                                             |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                                                                                                   |
| Cellular localization |                                                                                                                                                                                                                                                                                                                                                                   |
| Alternative Names     | DHRS1; Dehydrogenase/reductase SDR family member 1                                                                                                                                                                                                                                                                                                                |
| Background            | dehydrogenase/reductase 1(DHRS1) Homo sapiens This gene encodes a member of the short-chain dehydrogenases/reductases (SDR) family. The encoded enzyme contains a conserved catalytic domain and likely functions as an oxidoreductase. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Nov 2008], |