



## DHRS1 Polyclonal Antibody

Cat No: HR1AP5769

For research use only

### Overview

Product Name	DHRS1 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GenID?Human?	DHRS1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	DHRS1; Dehydrogenase/reductase SDR family member 1
Background	dehydrogenase/reductase 1(DHRS1) Homo sapiens This gene encodes a member of the short-chain dehydrogenases/reductases (SDR) family. The encoded enzyme contains a conserved catalytic domain and likely functions as an oxidoreductase. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Nov 2008],