

DGAT2L6 Polyclonal Antibody

Cat No: HR1AP6954

For research use only

Overview

Product Name	DGAT2L6 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	36kDa
GeneID?Human?	DGAT2L6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	DGAT2L6; DC3; Diacylglycerol O-acyltransferase 2-like protein 6; Diacylglycerol O-acyltransferase candidate 3; hDC3
Background	diacylglycerol O-acyltransferase 2 like 6(DGAT2L6) Homo sapiens This gene is a member of the diacylglycerol acyltransferase 2 family. The encoded protein is a putative acyltransferase and is most likely involved in the synthesis of di- or triacylglycerol, however its substrate specificity is currently unknown. [provided by RefSeq, Feb 2015],