



DHRS11 Polyclonal Antibody

Cat No: HR1AP8679

For research use only

Overview

Product Name	DHRS11 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GenID?Human?	DHRS11
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	DHRS11; Dehydrogenase/reductase SDR family member 11
Background	similarity:Belongs to the short-chain dehydrogenases/reductases (SDR) family.,subunit:Homotetramer.,

