



CYP4F11 Polyclonal Antibody

Cat No: HR1AP8048

For research use only

Overview

Product Name	CYP4F11 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	60kDa
GenelD?Human?	CYP4F11
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CYP4F11; Cytochrome P450 4F11; CYPIVF11
Background	cytochrome P450 family 4 subfamily F member 11(CYP4F11) Homo sapiens This gene, CYP4F11, encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This gene is part of a cluster of cytochrome P450 genes on chromosome 19. Another member of this family, CYP4F2, is approximately 16 kb away. Alternatively spliced transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008].