

COX19 Polyclonal Antibody

Cat No: HR1AP9057

For research use only

Overview

Product Name	COX19 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	COX19
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	COX19; Cytochrome c oxidase assembly protein COX19; hCOX19
Background	COX19, cytochrome c oxidase assembly factor(COX19) Homo sapiens COX19 encodes a cytochrome c oxidase (COX)-assembly protein. The S. cerevisiae Cox19 protein may play a role in metal transport to the mitochondrial intermembrane space and assembly of complex IV of the mitochondrial respiratory chain (Sacconi et al., 2005 [PubMed 16212937]).[supplied by OMIM, Mar 2008],