



CCS Polyclonal Antibody

Cat No: HR1AP9261

For research use only

Overview

Product Name	CCS Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	32kDa
GeneID?Human?	CCS
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CCS; Copper chaperone for superoxide dismutase; Superoxide dismutase copper chaperone
Background	copper chaperone for superoxide dismutase(CCS) Homo sapiens Copper chaperone for superoxide dismutase specifically delivers Cu to copper/zinc superoxide dismutase and may activate copper/zinc superoxide dismutase through direct insertion of the Cu cofactor. [provided by RefSeq, Jul 2008],