



V5-Tag Mouse mAb

Cat No: HR1AM2115

For research use only

Overview

Product Name	V5-Tag Mouse mAb
Source	Mouse
Applications	WB, IP, IF
Species Reactivity	N/A
Recommended Dilutions	WB 1:5,000 IP 1:200 IF 1:1,000
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	N/AkDa
GeneID?Human?	N/A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	N/A
Background	The V5 epitope tag is derived from a small epitope (Pk) present on the P and V proteins of the paramyxovirus of simian virus 5 (SV5). The V5 tag antibody can be helpful in detecting the recombinant proteins, some of which include transmembrane and secreted proteins fusion protein.