

AKR1CL1 Polyclonal Antibody

Cat No: HR1AP6931

For research use only

Overview

Product Name	AKR1CL1 Polyclonal Antibody
Source	Rabbit
Applications	IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	AKR1CL1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	AKR1CL1; Aldo-keto reductase family 1 member C-like protein 1
Background	similarity:Belongs to the aldo/keto reductase family.,